

Anti-Human/Mouse TRIB1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07241

3th Edition

Description: Interacts with Human, Mouse APK kinases and regulates activation of MAP kinases. May not display kinase activity.

Antigen: Synthetic peptide of human TRIB1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH 7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in most human tissues with the highest levels in skeletal muscle, thyroid gland, pancreas, peripheral blood leukocytes, and bone marrow.