

**Anti-Human/Mouse/Rat TM9SF1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA07246**

---

3th Edition

**Description:** Plays an essential role in autophagy.

**Antigen:** Synthetic peptide of human TM9SF1

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in lung, pancreas, kidney, liver, placenta, skeletal muscle, heart and brain. The amount in skeletal muscle, heart and brain were considerably lower than in the other tissues.