

Anti-Human RHCG Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07278

3th Edition

Description: Functions as an electroneutral and bidirectional ammonium transporter. May regulate transepithelial ammonia secretion.

Antigen: Synthetic peptide of human RHCG

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in brain, testis, placenta, pancreas, esophagus and prostate. Expressed in squamous epithelial tissues (at protein level). According to PubMed:11062476, specifically expressed in kidney.