

Anti-Human RHCE Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07280

3th Edition

Description:The Rh blood group system is the second most clinically significant of the blood groups, second only to ABO. It is also the most polymorphic of the blood groups, with variations due to deletions, gene conversions, and missense mutations. The Rh blood group includes this gene which encodes both the RhC and RhE antigens on a single polypeptide and a second gene which encodes the RhD protein. The classification of Rh-positive and Rh-negative individuals is determined by the presence or absence of the highly immunogenic RhD protein on the surface of erythrocytes.

Antigen:Synthetic peptide of human RHCE

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Restricted to tissues or cell lines expressing erythroid characters. Isoform 4g and isoform RhPI-Alpha are expressed in immature erythroblasts but not in mature erythroblasts.1 Publication