

**Anti-Human/Mouse PLXNB2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA07286**

---

3th Edition

**Description:**Members of the B class of plexins, such as PLXNB2 are transmembrane receptors that participate in axon guidance and cell migration in response to semaphorins (Perrot et al. (2002) [PubMed 12183458]).

**Antigen:**Synthetic peptide of human PLXNB2

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**