

**Anti-Human/Mouse GAS1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA07305**

---

3th Edition

**Description:** Growth arrest-specific 1 plays a role in growth suppression. GAS1 blocks entry to S phase and prevents cycling of normal and transformed cells. Gas1 is a putative tumor suppressor gene.

**Antigen:** Synthetic peptide of human GAS1

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**