

Anti-Human/Mouse/Rat CHGB Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07396

3th Edition

Description: This gene encodes a tyrosine-sulfated secretory protein abundant in peptidergic endocrine cells and neurons. This protein may serve as a precursor for regulatory peptides.

Antigen: Synthetic peptide of human CHGB

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in the adrenal medulla, and in pheochromocytoma. Not expressed in liver.