

Anti-Human/Mouse RAB17 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07411

3th Edition

Description:The Rab subfamily of small GTPases plays an important role in the regulation of membrane trafficking. RAB17 is an epithelial cell-specific GTPase.

Antigen:Synthetic peptide of human Rab17

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in melanocytes (at protein level).1 Publication