

Anti-Human/Mouse VPS37A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07419

3th Edition

Description: This gene belongs to the VPS37 family, and encodes a component of the ESCRT-I (endosomal sorting complex required for transport I) protein complex, required for the sorting of ubiquitinated transmembrane proteins into internal vesicles of multivesicular bodies. Expression of this gene is downregulated in hepatocellular carcinoma, and mutations in this gene are associated with autosomal recessive spastic paraplegia-53.

Antigen: Fusion protein of human VPS37A

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Widely expressed. Examined tissues include heart, brain, placenta, liver, skeletal muscle, kidney and pancreas. More abundant in liver. Strongly decreased or undetected in hepatomas.