

Anti-Human/Mouse USP39 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07466

3th Edition

Description: USP39 (ubiquitin specific peptidase 39), also known as SAD1, CGI-21, HSPC332 or SNRNP65, is a 565 amino acid nuclear protein belonging to the peptidase C19 family. Containing one UBP-type zinc finger, USP39 play a role in mRNA splicing and is part of a tri-snRNP complex. USP39 is thought to be required for the stabilization of the spindle checkpoint and may help support successful cytokinesis.

Antigen: Fusion protein of human USP39

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: