

**Anti-Human UGGT2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA07496**

---

3th Edition

**Description:**UDP-glucose:glycoprotein glucosyltransferase (UGT) is a soluble protein of the endoplasmic reticulum (ER) that selectively reglucosylates unfolded glycoproteins, thus providing quality control for protein transport out of the ER.

**Antigen:**Full length fusion protein

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Higher levels in kidney, pancreas, heart, and skeletal muscle.1 Publication