

Anti-Human UGGT2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07497

3th Edition

Description:UDP-glucose:glycoprotein glucosyltransferase (UGT) is a soluble protein of the endoplasmic reticulum (ER) that selectively reglucosylates unfolded glycoproteins, thus providing quality control for protein transport out of the ER.

Antigen:Full length fusion protein

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Higher levels in kidney, pancreas, heart, and skeletal muscle.1 Publication