

**Anti-Human/Mouse/Rat UBAP1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA07509**

---

3th Edition

**Description:** This gene is a member of the UBA domain family, whose members include proteins having connections to ubiquitin and the ubiquitination pathway. The ubiquitin associated domain is thought to be a non-covalent ubiquitin binding domain consisting of a compact three helix bundle. This particular protein originates from a gene locus in a refined region on chromosome 9 undergoing loss of heterozygosity in nasopharyngeal carcinoma (NPC).

**Antigen:** Fusion protein of human UBAP1

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Ubiquitous. Highly expressed in heart, brain, placenta, lung, liver, skeletal muscle and pancreas.1  
Publication