

Anti-Human/Mouse TXNDC5 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07513

3th Edition

Description: This locus represents naturally occurring read-through transcription between the neighboring MUTED (muted homolog) and TXNDC5 (thioredoxin domain containing 5) genes on chromosome 6. The read-through transcript is a candidate for nonsense-mediated mRNA decay (NMD) and is unlikely to produce a protein product.

Antigen: Full length fusion protein

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: