

Anti-Human/Mouse TUSC2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07520

3th Edition

Description: This gene is a highly conserved lung cancer candidate gene. No other information about this gene is currently available.

Antigen: Full length fusion protein

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Strong expression in heart, lung, skeletal muscle, kidney, and pancreas, followed by brain and liver, lowest levels in placenta.