

Anti-Human TULP2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07521

3th Edition

Description: TULP2 is a member of a family of tubby-like genes (TULPs) that encode proteins of unknown function. Members of this family have been identified in plants, vertebrates, and invertebrates. The TULP proteins share a conserved C-terminal region of approximately 200 amino acid residues.

Antigen: Fusion protein of human TULP2

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Strongly expressed in testis. Also expressed in retina. Expressed in cancer cell lines.