

**Anti-Human TULP2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA07522**

---

3th Edition

**Description:** TULP2 is a member of a family of tubby-like genes (TULPs) that encode proteins of unknown function. Members of this family have been identified in plants, vertebrates, and invertebrates. The TULP proteins share a conserved C-terminal region of approximately 200 amino acid residues.

**Antigen:** Fusion protein of human TULP2

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Strongly expressed in testis. Also expressed in retina. Expressed in cancer cell lines.