

Anti-Human/Mouse/Rat TUBA1C Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07523

3th Edition

Description: Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.

Antigen: Fusion protein of human TUBA1C

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: