



Anti-Human TRAK1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07558

3th Edition

Description:Involved in the regulation of endosome-to-lysosome trafficking, including endocytic trafficking of EGF-EGFR complexes and GABA-A receptors.

Antigen:Fusion protein of human TRAK1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:High expression in spinal cord and moderate expression in all other tissues and specific brain regions examined. Expressed in all cell lines examined.1 Publication