

Anti-Human Tmprss11a Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07571

3th Edition

Description: Probable serine protease which may play a role in cellular senescence. Overexpression inhibits cell growth and induce G1 cell cycle arrest.

Antigen: Fusion protein of human Tmprss11a

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in esophagus, liver, colon and lung. Down-regulated in esophagus cancers.