

Anti-Human TMPRSS11E Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07572

3th Edition

Description: Serine protease which possesses both gelatinolytic and caseinolytic activities. Shows a preference for Arg in the P1 position (By similarity).

Antigen: Fusion protein of human TMPRSS11E

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expression can only be detected in tissues derived from the head and neck, and in skin, prostate and testis.