

Anti-Human/Mouse/Rat THAP4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07619

3th Edition

Description: THAP domain-containing protein 4 is a protein that in humans is encoded by the THAP4 gene.

Antigen: Fusion protein of human THAP4

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: