

Anti-Human/Mouse/Rat STX12 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07640

3th Edition

Description: Syntaxin 12, also known as STX12, STX13 or STX14, is a 276 amino acid single-pass membrane protein that contains one t-SNARE coiled-coil homology domain and belongs to the Syntaxin family. Syntaxin 12 regulates protein transport between late endosomes and the trans-Golgi network and interacts with ABC1 (ATP-binding cassette transporter A1), a protein that facilitates cellular release of choline-phospholipids and cholesterol to apoA-I (apolipoprotein A-I).

Antigen: Full length fusion protein

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: