

Anti-Human/Mouse SRPK3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07641

3th Edition

Description: This gene encodes a protein kinase similar to a protein kinase which is specific for the SR (serine/arginine-rich domain) family of splicing factors. A highly similar protein has been shown to play a role in muscle development in mice. Multiple transcript variants encoding different isoforms have been found for this gene.

Antigen: Fusion protein of human SRPK3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Exclusively expressed in skeletal and heart muscle.