

Anti-Human SERPINA12 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07659

3th Edition

Description: Adipokine that modulates insulin action by specifically inhibiting its target protease KLK7 in white adipose tissues.

Antigen: Fusion protein of human SERPINA12

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in visceral adipose tissues.