

Anti-Human/Mouse SEMA4C Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07666

3th Edition

Description: Cell surface receptor for PLXNB2 that plays an important role in cell-cell signaling. PLXNB2 binding promotes downstream activation of RHOA and phosphorylation of ERBB2 at 'Tyr-1248'. Required for normal brain development, axon guidance and cell migration (By similarity). Probable signaling receptor which Human, Mouse play a role in myogenic differentiation through activation of the stress-activated MAPK cascade.

Antigen: Fusion protein of human SEMA4C

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: