

**Anti-Human/Mouse/Rat SDF4 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA07673**

---

3th Edition

**Description:** This gene encodes a stromal cell derived factor that is a member of the CREC protein family. The encoded protein contains six EF-hand motifs and calcium-binding motifs. This protein localizes to the Golgi lumen and may be involved in regulating calcium dependent cellular activities.

**Antigen:** Fusion protein of human SDF4

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Ubiquitous. Isoform 5 is expressed in pancreas.