

Anti-Human TRPM6 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07733

3th Edition

Description: This gene is predominantly expressed in the kidney and colon, and encodes a protein containing an ion channel domain and a protein kinase domain. It is crucial for magnesium homeostasis, and plays an essential role in epithelial magnesium transport and in the active magnesium absorption in the gut and kidney. Mutations in this gene are associated with hypomagnesemia with secondary hypocalcemia. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

Antigen: Synthetic peptide of human TRPM6

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Highly expressed in kidney and colon. Isoform TRPM6a and isoform TRPM6b, are coexpressed with TRPM7 in kidney, and testis, and are also found in several cell lines of lung origin. Isoform TRPM6c is detected only in testis and in NCI-H510A small cell lung carcinoma cells.