

Anti-Human/Mouse/Rat PFN4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07763

3th Edition

Description: This gene belongs to the profilin family. It binds to actin and affects the structure of the cytoskeleton. Mostly expression in testes, pituitary gland, midbrain, bone marrow, vulva.

Antigen: Recombinant protein of human PFN4

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Detected in testis (at protein level).