

Anti-Human UTS2R Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07826

3th Edition

Description: This gene codes high affinity receptor for urotensin-2 and urotensin-2B. The activity of this receptor is mediated by a G-protein that activate a phosphatidylinositol-calcium second messenger system

Antigen: Synthetic peptide of human UTS2R

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Most abundant expression in the heart and pancreas.