

Anti-Human/Mouse EVA1A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07863

3th Edition

Description: TMEM166, also known as FAM176A (family with sequence similarity 176, member A), is a 152 amino acid protein encoded by a gene mapping to human chromosome 2. The second largest human chromosome, 2 consists of 237 million bases encoding over 1,400 genes and making up approximately 8% of the human genome. A number of genetic diseases are linked to genes on chromosome 2. Harlequin ichthyosis, a rare and morbid skin deformity, is associated with Human, Mouse mutations in the ABCA12 gene.

Antigen: Synthetic peptide of human EVA1A

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in lung, kidney, liver, pancreas, placenta, but not in heart and skeletal muscle.