

**Anti-Human/Mouse/Rat TRAF3IP1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA07906**

---

3th Edition

**Description:** This protein plays an inhibitory role on IL13 signaling by binding to IL13RA1. Involved in suppression of IL13-induced STAT6 phosphorylation, transcriptional activity and DNA-binding. Recruits TRAF3 and DISC1 to the microtubules. It can interact with IL13RA1 and binds to microtubules, TRAF3 and DISC1.

**Antigen:** Recombinant protein of human TRAF3IP1

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Ubiquitous.