

Anti-Human/Mouse KLF5 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07919

3th Edition

Description: This gene encodes a member of the Kruppel-like factor subfamily of zinc finger proteins. Since the protein localizes to the nucleus and binds the epidermal growth factor response element, it is thought to be a transcription factor.

Antigen: Synthetic peptide of human KLF5

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed only in testis and placenta.