

Anti-Human/Mouse BPIFB2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07964

3th Edition

Description: This gene encodes a member of the lipid transfer/lipoplysaccharide binding protein (LT/LBP) gene family. It is highly expressed in hypertrophic tonsils. This gene and three other members of the LT/LBP gene family form a cluster on the long arm of chromosome 20.

Antigen: Synthetic peptide of human BPIFB2

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Highly expressed in tonsils, especially in hypertrophic tonsils. Detected at very low levels in fetal liver.