

**Anti-Human/Mouse STK32C Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA07968**

---

3th Edition

**Description:** The protein encoded by this gene is a member of the serine/threonine protein kinase family. The specific function of this kinase is not known.

**Antigen:** Synthetic peptide of human STK32C

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**