

Anti-Human/Mouse STAB1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07980

3th Edition

Description: This gene encodes a large, transmembrane receptor protein which Human, Mouse may function in angiogenesis, lymphocyte homing, cell adhesion, or receptor scavenging. The protein contains 7 fasciclin, 16 epidermal growth factor (EGF)-like, and 2 laminin-type EGF-like domains as well as a C-type lectin-like hyaluronan-binding Link module. The protein is primarily expressed on sinusoidal endothelial cells of liver, spleen, and lymph node.

Antigen: Synthetic peptide of human STAB1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: High levels found in spleen, lymph node, liver and placenta. Also expressed in endothelial cells.