

**Anti-Human/Mouse/Rat SRGAP2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA07988**

---

3th Edition

**Description:** This locus encodes a member of the SLIT-ROBO Rho GTPase activating protein family. The encoded protein stimulates GTPase activity of Rac1, and plays a role in cortical neuron development. This locus has several paralogs on human chromosome 1 resulting from segmental duplication. While this locus itself is conserved among various species, the paralogs are found only in the genus Homo, and not in the genomes of non-human great apes.

**Antigen:** Synthetic peptide of human SRGAP2

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expression at low levels in placenta, kidney, and ovary, and little or no expression in all other tissues tested.