

Anti-Human A4GNT Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07996

3th Edition

Description: This gene encodes a protein from the glycosyltransferase 32 family. The enzyme catalyzes the transfer of N-acetylglucosamine (GlcNAc) to core 2 branched O-glycans. It forms a unique glycan, GlcNAc α 1 \rightarrow 4Gal β 1 \rightarrow 6R and is largely associated with the Golgi apparatus membrane.

Antigen: Recombinant protein of human A4GNT

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1 μ g/ml.

Specificity: Detected in stomach and pancreas.