

**Anti-Human/Mouse/Rat AASDHPPT Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA08005**

---

3th Edition

**Description:**The protein encoded by this gene is similar to *Saccharomyces cerevisiae* LYS5, which is required for the activation of the alpha-aminoadipate dehydrogenase in the biosynthetic pathway of lysine. Yeast alpha-aminoadipate dehydrogenase converts alpha-biosynthetic-aminoadipate semialdehyde to alpha-aminoadipate. It has been suggested that defects in the human gene result in pipecolic acidemia.

**Antigen:**Recombinant protein of human AASDHPPT

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Detected in heart, skeletal muscle, placenta, testis, brain, pancreas, liver and kidney.