

Anti-Human/Mouse/Rat SNX3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08036

3th Edition

Description: This gene encodes a member of the sorting nexin family. Members of this family contain a phosphoinositide binding domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region, like most family members. This protein interacts with phosphatidylinositol-3-phosphate, and is involved in protein trafficking. A pseudogene of this gene is present on the sex chromosomes. Alternative splicing results in multiple transcript variants encoding distinct isoforms.

Antigen: Synthetic peptide of human SNX3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: