

Anti-Human ACOT1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08057

3th Edition

Description: Acyl-CoA thioesterases such as ACOT1, hydrolyze acyl-CoAs to the free fatty acid and CoA. ACOTs therefore play key roles in maintaining the intracellular ratio between CoA esters of various lipids and free fatty acids. Acyl-CoA thioesterases are a group of enzymes that catalyze the hydrolysis of acyl-CoAs to the free fatty acid and coenzyme A (CoASH), providing the potential to regulate intracellular levels of acyl-CoAs, free fatty acids and CoASH. Active towards fatty acyl-CoA with chain-lengths of C12-C16.

Antigen: Recombinant protein of human ACOT1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: