

Anti-Human/Mouse/Rat SLIT3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08070

3th Edition

Description:The protein encoded by this gene is secreted, likely interacting with roundabout homolog receptors to effect cell migration. Two transcript variants encoding different isoforms have been found for this gene.

Antigen:Synthetic peptide of human SLIT3

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Predominantly expressed in thyroid.