

Anti-Human/Mouse/Rat MYL9 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08085

3th Edition

Description: Myosin, a structural component of muscle, consists of two heavy chains and four light chains. The protein encoded by this gene is a myosin light chain that may regulate muscle contraction by modulating the ATPase activity of myosin heads. The encoded protein binds calcium and is activated by myosin light chain kinase. Two transcript variants encoding different isoforms have been found for this gene.

Antigen: Synthetic peptide of human MYL9

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Smooth muscle tissues and in some, but not all, nonmuscle cells.