

**Anti-Human SLC25A38 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA08139**

---

3th Edition

**Description:** This gene is a member of the mitochondrial carrier family. The encoded protein is required during erythropoiesis and is important for the biosynthesis of heme. Mutations in this gene are the cause of autosomal congenital sideroblastic anemia.

**Antigen:** Synthetic peptide of human SLC25A38

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Preferentially expressed in erythroid cells.