

Anti-Human/Mouse/Rat SLC8B1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08142

3th Edition

Description:SLC24A6 belongs to a family of potassium-dependent sodium/calcium exchangers that maintain cellular calcium homeostasis through the electrogenic countertransport of 4 sodium ions for 1 calcium ion and 1 potassium ion (Cai and Lytton, 2004 [PubMed 14625281]).

Antigen:Synthetic peptide of human SLC24A6

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Present in pancreatic beta-cells (at protein level).1 Publication