

**Anti-Human/Mouse ALG9 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA08222**

---

3th Edition

**Description:** This gene encodes an alpha-1,2-mannosyltransferase enzyme that functions in lipid-linked oligosaccharide assembly. Mutations in this gene result in congenital disorder of glycosylation type II. Multiple transcript variants encoding different isoforms have been found for this gene.

**Antigen:** Recombinant protein of human ALG9

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Ubiquitously expressed; with highest levels in heart, liver and pancreas.