

Anti-Human/Mouse/Rat CWC27 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08237

3th Edition

Description:CWC27 belongs to the cyclophilin-type PPlase family and contains 1 PPlase cyclophilin-type domain. There are two isoforms. PPlases accelerate the folding of proteins.

Antigen:Synthetic peptide of human CWC27

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: