

Anti-Human/Mouse/Rat RNF126 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08248

3th Edition

Description:RNF126 contains a RING type zinc finger domain, a motif known to be involved in protein-protein and protein-DNA interactions. The specific function of this protein is unknown. There are two named isoforms.

Antigen:Synthetic peptide of human RNF126

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Highly expressed in liver and testis.