

Anti-Human SLC12A6 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08269

3th Edition

Description: This gene is a member of the K-Cl cotransporter (KCC) family. K-Cl cotransporters are integral membrane proteins that lower intracellular chloride concentrations below the electrochemical equilibrium potential. The proteins encoded by this gene are activated by cell swelling induced by hypotonic conditions. Alternate splicing results in multiple transcript variants encoding different isoforms. Mutations in this gene are associated with agenesis of the corpus callosum with peripheral neuropathy.

Antigen: Synthetic peptide of human SLC12A6

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Highly expressed in heart, brain and kidney. Detected at lower levels in skeletal muscle, placenta, lung and pancreas. Detected in umbilical vein endothelial cells. Isoform 2 is more abundant in kidney. Isoform 5 is testis specific. Expressed in the proximal tubule of the kidney (at protein level).