

Anti-Human MAP4K1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA08308

3th Edition

Description: This gene belongs to the protein kinase superfamily, STE Ser/Thr protein kinase family, STE20 subfamily. It plays a role in the response to environmental stress and hematopoietic lineage decisions and growth regulation. Appears to act upstream of the JUN N-terminal pathway.

Antigen: Synthetic peptide of human MAP4K1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed primarily in hematopoietic organs, including bone marrow, spleen and thymus. Also expressed at very low levels in lung, kidney, mammary glands and small intestine.